



Laboratório de Genómica Funcional de Plantas

Docente: Andreia Figueiredo

cDNA synthesis protocol

Important: For qPCR you have to be sure that total RNA is not contaminated with gDNA.

Normally 2-5ug of total RNA is used for cDNA synthesis.

In grapevine normally we use 2.5ug of tRNA.

Reagents:

- Oligo dt primer (18 or 23 oligomer) / hexamers or decamers
- Buffer RT 5x
- Reverse transcriptase 200U/ul
- DNTPs (10mM)
- Ribolock RNase inhibitor 40 U/ul

Equipment:

- PCR

Protocol:

- total reaction volume: 20ul
- Add 1 ul of oligo dt primer to 2.5 ug of tRNA in a PCR tube and incubate 10min at 70°C on a thermocycler.
- put the PCR tube on ice
- Add to the PCR tube the following:
 - 5x RT buffer ----- 4 ul
 - 10mM DNTPs ----- 2 ul
 - Ribolock ----- 0.5 ul
 - Reverse Transcriptase -- 0.5 ul



-Optional: you may incubate the mix at 37 °C during 5 min on a thermocycler prior to the addition of the Reverse Transcriptase

-Incubate at 42 °C for 60 min (up to 90 min)

-Stop the reaction by incubating at 70 °C for 10 min

-Keep the cDNA at -20 °C until further use

Check for cDNA on agarose gel

Run 1ul of the synthesized cDNA on 1-1.2% agarose gel. The cDNA should look like a smear on gel.