




Ecotoxicology TP Course

Concepts, Tests & Biomarkers

Bernardo Duarte

Faculty of Sciences of the University of Lisbon
MARE – Marine and Environmental Sciences Centre



TP7 Effect and defense biomarkers



Control



10 ug/L Glifosato



250 ug/L Glifosato




500 ug/L Glifosato


- Add 1.5 mL TCA 20% (w/v) to the pellet.
- Heat the sample for 20 min at 100 °C.
- Take 1 mL of the TCA sample and add it to 1 mL TBA 0.5 %.
- Heat the sample for an additional 20 min at 100 °C.
- Read the absorbance at 532 nm e 600 nm (ϵ [MDA] = 155 mM⁻¹ cm⁻¹):

$$A_{532 \text{ nm}} - A_{600 \text{ nm}} = [\text{MDA}]_{\text{mM}} \times \epsilon_{\text{MDA}}$$

[Glyphosate]	[TBARS] pmol cell 10 ⁻⁶
0	18,1364
10	35,18561
250	127,7721
500	158,1783





[Glyphosate]	[TBARS] pmol cell 10 ⁻⁶	Enhancement (%)
0	18,1364	
10	35,18561	
250	127,7721	
500	158,1783	





[Glyphosate]	[TBARS] pmol cell 10 ⁻⁶	Enhancement (%)
0	18,1364	0%
10	35,18561	94%
250	127,7721	605%
500	158,1783	772%

EC₅₀ = ?



[Glyphosate]	[TBARS] pmol cell 10 ⁻⁶	Enhancement (%)
0	18,1364	0%
10	35,18561	94%
250	127,7721	605%
500	158,1783	772%



EC₅₀ = 6.22 mg/L



OXIDATIVE STRESS:

- Oxidative stress reflects an **imbalance** between the systemic **manifestation of reactive oxygen species** and a **biological system's ability to readily detoxify** the reactive intermediates or to repair the resulting damage.
- Disturbances in the normal redox state of cells can cause toxic effects through the production of peroxides and free radicals that damage all components of the cell, including proteins, lipids, and DNA.
- Oxidative stress from oxidative metabolism causes base damage, as well as strand breaks in DNA. Base damage is mostly indirect and caused by reactive oxygen species (ROS) generated, e.g. O_2^- (superoxide radical), OH (hydroxyl radical) and H_2O_2 (hydrogen peroxide).
- Further, some reactive oxidative species act as cellular messengers in redox signalling. Thus, oxidative stress can cause disruptions in normal mechanisms of cellular signalling.
- **IMPORTANT:** ROS are generated during the normal cell functioning thus the cells also developed anti-oxidant mechanisms to counteract the negative effects produced by ROS.
- In Ecotoxicology, the production or activity of these ROS counteractive measures as well as the ROS-biomolecules reaction products are used to evaluate the cell oxidative stress level.

SOURCES OF OXIDATIVE STRESS



PHASE I: BIOTRANSFORMATION


Phase I enzymes and mechanisms transform xenobiotics into less harmful molecules, but that may have a ROS-generating potential (for e.g.: CYP450 enzymes).

PHASE II: CONJUGATION


Phase II mechanisms are composed by enzymatic and non-enzymatic **anti-oxidant mechanisms** that work in conjugation to quench the ROS generated directly by the xenobiotic or by the Phase I biotransformed xenobiotic.

PHASE III: EXCRETION


Phase III mechanisms are based in **membrane proteins that can excrete** directly or throughout vesicle compartments the transformed and/or inactivated xenobiotic to the extracellular environment.



Oxygen	O_2
Superoxide anion	O_2^-
Peroxide anion	O_2^{2-}
Hydrogen peroxide	H_2O_2
Hydroxyl radical	OH
Hydroxyl anion	OH^-
Singlet Oxygen	1O_2
Hypochloric acid	$HOCl$



FREE RADICAL TOXICITY



DNA AND PROTEIN DAMAGE

ROS covalent binding to DNA leads to double strand disruption.
ROS interaction with Proteins induces protein oxidation/carboxylation.

ENZYMATIC DEFENSE SYSTEM

Superoxide dismutase isoforms and several peroxidases decompose ROS into harmless substances.

NON-ENZYMATIC ANTIOXIDANT

Vitamins, phenolics, flavonoids, thiol molecules and other anti-oxidant are able to quench directly ROS molecules into stable harmless radicals.


MEMBRANE DAMAGE

ROS interaction membrane fatty acids induces the formation of lipid hydroperoxides and lipid radicals inducing membrane disruption.


BIOMARKER


A trait or molecular entity that can be measured experimentally and indicate the occurrence of a certain function (normal or pathological) of a certain organisms towards a specific stressor.

Ideally these biomarkers should produce a dose related response towards the stressor applied.



BIOMARKERS : ENZYME EXTRACTION





- All these enzymes are part of the normal functioning of the cells.
- Cells produce ROS as part of their normal cellular metabolism.
- In Ecotoxicology, the activity of these enzymes towards the basal cell functioning is evaluated under stress conditions as a measure of the cell stress.

ASCORBATE PEROXIDASE (Apx, *ONLY PLANTS*)

Uses ascorbate molecules to quench the oxidative power from H_2O_2 , converting hydrogen peroxide into H_2O .

GLUTATHIONE PEROXIDASE (Gpx) AND REDUCTASE (GR)

Glutathione peroxidase used reduced glutathione molecules to convert H_2O_2 into H_2O . The oxidized form of glutathione (GSSH) is then reduced back by glutathione reductase at the expense of NADH.

CATALASE (CAT)

Converts hydrogen peroxide into water.

PEROXIDASE (POX)

A class of several peroxidasic enzymes that convert H_2O_2 into H_2O or OH^\cdot at the expense of NAD(P)H.

SUPEROXIDE DISMUTASE (SOD)

This enzyme converts the superoxide anion into hydrogen peroxide feeding the peroxidasic system.

BIOMARKERS : SUPEROXIDE DISMUTASE ASSAY (PHASE II)

PRINCIPLES OF THE METHOD:

- Pyrogallol (Pyr) has a absorption peak at 325 nm.
- Pyrogallol has an auto-oxidation (degradation) rate.
- The method measures the inhibition of oxidized Pyrogallol for 120 seconds at 325 nm by Superoxide Dismutase.



BLANK CELL


1. 100 µg protein (calculate the extract volume needed).
2. Remaining volume with PBS buffer.

To zero the spectrophotomer.



SOD ASSAY

1. 100 µg protein (calculate the extract volume needed).
2. Remaining volume with PBS buffer.
3. 0.3 mM Pyrogallol



CALCULATIONS

$$\text{Slope}_{\text{BLANK}} = (\text{Abs}_{t=120} - \text{Abs}_0)_{\text{BLANK}} / (120 - 0)$$

$$\text{Slope}_{\text{SOD}} = (\text{Abs}_{t=120} - \text{Abs}_0)_{\text{BLANK}} / (120 - 0)$$

$$\% \text{ INHIBITION PYROGALLOL} = \text{Slope}_{\text{SOD}} / \text{Slope}_{\text{BLANK}}$$

$$\text{SOD (U/mL)} = \% \text{ INHIBITION PYROGALLOL} \times 50\%$$

BIOMARKERS : SUPEROXIDE DISMUTASE ASSAY (PHASE II)



1. Extract the proteins with Extraction PBS Buffer (0.5 mL)
2. Sonicate
3. Centrifuge to precipitate cellular debris
4. Use the extract for spectrophotometric assay (60 sec at 325 nm)

	Blank	Reaction
Ultra-pure Water	360 μL	360 μL
Buffer	630 μL	550 μL
Extract	10 μL	10 μL
Pyrogallol	0 μL	80 μL