

Bioinorganic Coordination Chemistry: Copper(II) Tetraphenylporphyrinate

Note: This experiment requires 2 hours for Part A and 6 hours for Part B.

Metal ions play vital roles in many biological processes, and at least seven transition metals (iron, zinc, copper, manganese, cobalt, nickel, and molybdenum) are essential to almost all life on earth. These metals are key components of many important proteins. In some cases, the metals coordinate to the nitrogen, sulfur, or oxygen atoms in the side chains of certain amino acids that make up the protein's structure; among these "metal-binding" amino acids are histidine, cysteine, methionine, tyrosine, aspartic acid, and glutamic acid. In other cases, however, the transition metals are bound to special ligands, the most important of which are the porphyrins. Metal-bearing porphyrin complexes are called metalloporphyrins.

Porphyrins have the general structure shown in Figure 23-1; they are compounds with a central 16-membered ring consisting of four pyrrole subunits linked by one-carbon bridges. The porphyrin ring is polyunsaturated and completely conjugated; consequently, porphyrins and their complexes with transition metals are intensely colored. In metalloporphyrins, a metal atom coordinates




Figure 23-1

General structure of a porphyrin with peripheral groups R₁–R₁₂; structure of the copper(II) complex Cu(TPP).


to the four nitrogen atoms and displaces the two central hydrogen atoms. One important metalloporphyrin complex that is found in almost all animals is hemoglobin, which contains four iron-porphyrin units. In vertebrates, hemoglobin is responsible for the transport of O_2 from the lung to cells throughout the body. Metalloporphyrin complexes perform a variety of other important biochemical functions; for example, they serve as electron-transfer relays and as oxidation catalysts.

Closely related to metalloporphyrins are the chlorins, corrins, and corphins, which have similar though not identical structures. Two important examples of such porphyrin-like molecules are chlorophyll and vitamin B_{12} (see Experiment 21). Phthalocyanines, which are nitrogen-rich analogues of porphyrins, are produced and used industrially as pigments and catalysts.

Naturally occurring porphyrins generally have a variety of different organic groups on the periphery (exterior) of the ring. The synthesis of such polysubstituted porphyrin rings is a challenging task for the chemist. For many purposes, however, the characteristic chemical properties associated with metalloporphyrins are exhibited by simpler analogues with small peripheral groups. The most important examples of such synthetic analogues are the complexes of *meso*-tetraphenylporphyrin (abbreviated H_2TPP), where the "meso" designation means that the phenyl groups are located on the four carbon atoms that bridge between the pyrrole rings. All the other peripheral groups in H_2TPP are hydrogen atoms, and the chemical formula of H_2TPP is $C_{44}H_{30}N_4$. In this experiment, you will prepare H_2TPP and convert it to its copper complex $Cu(TPP)$, whose structure is shown in Figure 23-1.

This experiment begins with the preparation of H_2TPP by the condensation of four molecules each of benzaldehyde and pyrrole. This reaction does not give a high yield but the starting materials are inexpensive and the product is easily isolated. The low yields of this reaction illustrate the difficulty of assembling a large ring in a "one-pot" reaction. In some cases, but not this one, the yields of such cyclization reactions can be improved by adding metal ions to the solution, which assist in the assembly of the macrocycle by binding to the reactants and orienting them in such a fashion to favor formation of the ring.

The stoichiometry of the reaction of benzaldehyde and pyrrole is as follows:



The mechanism of the cyclization reaction is known in some detail: It involves formation of a carbocation by addition of a proton to benzaldehyde, followed by electrophilic attack of the carbocation at the α position of pyrrole. Loss of water generates a new carbocation, which attacks a second pyrrole ring. These steps are repeated, and eventually a nonconjugated macrocycle is formed. Oxidation of this macrocycle by O_2 generates the fully conjugated porphyrin ring.

Because the oxidation step is not quantitative, however, one of the principal contaminants is *meso*-tetraphenylchlorin (H_2TPC), a compound that contains two more hydrogen atoms than H_2TPP . The chemical structure of H_2TPC is shown in Figure 23-2; the extra hydrogen atoms are present on the β -carbon atoms of one of the pyrrole rings:





Figure 23-2
Structure of *meso*-tetraphenylchlorin (H_2TPC).

Tetraphenylporphyrin and its metal complexes are best separated from tetraphenylchlorin impurities by chromatography.

The conversion of H_2TPP to $Cu(TPP)$ is achieved by the addition of an excess of copper(II) acetate. The copper(II) center replaces the two nitrogen-bound protons at the center of the H_2TPP ring:



The insertion of a metal ion into the porphyrin ring is sometimes very slow even though the reaction is very favorable thermodynamically. This somewhat surprising behavior reflects the rigidity and steric constraints of the macrocycle, which make it difficult for the porphyrin nitrogen atoms to approach a metal center that is already surrounded by other ligands (such as water or acetate groups).

Thin-Layer and Column Chromatography

One of the goals of this experiment is to illustrate standard chromatographic techniques that are used in the isolation of pure compounds from complex reaction mixtures. Thin-layer chromatography (TLC) will be used to test the purity of small amounts of the $Cu(TPP)$ sample and to explore the chromatographic conditions necessary to purify it. Once these conditions are established, column chromatography will be utilized to separate larger amounts of $Cu(TPP)$ from copper(II) tetraphenylchlorin and other impurities.

Separations effected by both TLC and column chromatography are based on the tendency of molecules to bind to certain solids called adsorbents. The adsorbents most frequently used are silica gel and alumina. Silica gel, which will be used as the adsorbent in this experiment, is to a first approximation a hydrated form of silicon dioxide, $SiO_2 \cdot x H_2O$. It does, however, contain significant amounts of other inorganic salts whose amounts vary from one silica gel preparation to another. For TLC, the silica gel is often mixed with plaster to help bind the gel to a glass support.

When silica gel is heated in strongly acidic or basic solutions, it acquires acidic or basic properties. Acid-treated silica gel strongly adsorbs (or binds) basic compounds such as amines, whereas base-treated silica gel adsorbs acidic compounds. The adsorption properties of silica gels also depend on their water content. If the gels are strongly heated under vacuum, water is driven off the silica gel leaving sites where other polar molecules strongly adsorb. Less strongly adsorbing silica gel can be prepared by adding back small amounts of H_2O to occupy some of the adsorption sites. By altering the water content, it is therefore possible to control the degree to which silica gel binds various compounds.

In TLC, a thin layer (0.1–2 mm thick) of the adsorbent is spread onto a flat surface. The TLC plates can be purchased commercially and can also be prepared by coating microscope slides. A small amount of the sample to be separated is dissolved in a small volume of a suitable solvent. It is essential that the sample to be tested be completely dissolved. With a capillary, a spot (3–5 mm in diameter) of the solution is placed about 8 mm from one end of the TLC plate. The plate is allowed to dry and then the procedure is repeated to add sample to the same spot. The spot should be kept small for maximum separation of the components. The TLC plate is then placed in a bottle that contains a few milliliters of solvent (Fig. 23-3). The solvent level must be lower than the spot on the




Figure 23-3
Separation using TLC.

plate. The bottle is capped and left undisturbed while the solvent rises up the silica gel by capillary action. If the proper solvent is chosen, the mixture will begin to move up the plate behind the solvent front. Ideally, different compounds present in the mixture will move up the plate at different rates, and thus will be separated.

When the solvent has moved about three-fourths of the way up, the plate should be removed from the bottle. If the compounds in the mixture are colored, it will be obvious if a separation has occurred. If some of the compounds are colorless, their locations on the TLC plate can usually be established by placing the air-dried TLC plate in a bottle containing a few crystals of iodine. The iodine sublimates and adsorbs in the areas where the compounds are located. Thus, dark brown spots on the plate indicate the locations of the components of the original sample. Alternatively, some TLC plates contain a fluorescent indicator:

Upon illumination with an ultraviolet (UV) light, the air-dried plate will glow in all places except the spots where compounds are located.

The separation of mixtures into their components by chromatography mainly depends on differences in the adsorption tendencies and solubilities of the components. Compounds that are weakly adsorbed and that are readily soluble in the solvent will move or elute quickly; in contrast, compounds that are strongly adsorbed and that are poorly soluble in the solvent will elute slowly. Finding the solvent that effects the best separation of the components is not easy, but generally a solvent that dissolves the desired compound moderately well will allow the compound to move up the plate; it can only be hoped that the impurities do not migrate at the same rate as does the compound of interest. If the solvent that was chosen does not separate the components of the mixture, other solvents either more or less polar than the first should be tried until a solvent that gives a separation is found. The polarities of some common chromatographic solvents increase in the following order:

↑
increasing polarity
↓

alkanes
toluene
benzene
dichloromethane
chloroform
diethyl ether
ethyl acetate
acetone
ethanol
methanol
water

It is sometimes convenient to use mixtures of solvents. For example, mixtures of ethyl acetate in dichloromethane often succeed in effecting a useful separation where other solvents do not.

At this point it is probably obvious that the successful choice of adsorbent and solvent is an art that is learned largely by doing chromatographic separations. The references at the end of this experiment do offer, however, many hints on how to use these techniques effectively.

Having established the solvent or solvent mixture that will separate the sample on TLC plates, it is hoped that the same solvent can be used to separate larger quantities of the sample on a silica gel chromatography column. Generally, this is possible. It is necessary, however, to use a much larger silica gel particle size for column chromatography (80–200 mesh) than that used in TLC (finer than 200 mesh).

The chromatography column (Fig. 23-4) consists of a glass tube that typically is about 3 cm in diameter and 30 cm long. First, a small glass wool plug is pushed to the bottom of the column, and a 5-mm layer of sand is added. A slurry of the silica gel in the solvent to be used in the separation is then poured onto the sand. The column is drained until the solvent level is the same as the top of the silica gel, and then a mixture of the sample and silica gel (both suspended in a few milliliters of the solvent) is added to the column. (Alternatively, the sample may




Figure 23-4
Chromatography column.

be dissolved in a small volume of the solvent and added to the column.) A 5-mm layer of sand is then added to the top of the silica gel.

From this point on, the solvent level should never fall below the top of the silica gel, because channels in the column will result, and the solution will pass down the channels without properly percolating through the adsorbent. The solvent is added to the column, and the movement of the compounds down the column begins. This process is called *elution*. The rate at which solvent is passed through the column is called the *flow-rate*; slow flow-rates give better separations than do fast flow-rates.

If the compounds in the sample are colored, it is easy to determine which fractions should be collected to obtain the desired products. If the compounds are colorless, fractions may have to be collected at regular volume intervals and examined for the presence of the desired compounds spectrophotometrically or by other techniques. Whereas the initial eluting solvent may elute one or more of the compounds, other compounds may require (as noted by TLC) more polar solvents to move them down the column. If a change in solvent is required, it is best to introduce it gradually by using first a mixture of the initial and the subsequent solvent and then finally the pure new solvent. Sometimes an abrupt change of solvent leads to the evolution of large amounts of heat when the new solvent adsorbs to the silica gel. Thermal expansion of the solvent creates channels in the adsorbent that can destroy the efficiency of the column.

The fractions eluted from the column that contain the desired compounds may simply be evaporated to dryness to give the pure product. Evaporation sometimes does not give a crystalline solid, and in such cases recrystallization of the material usually gives purer and better looking product.

Column chromatography has numerous variations. Although silica gel and alumina (Al_2O_3) are the most common adsorbents, many others have also been used. For materials that decompose at room temperature, chromatographic separations can be carried out in cooled, jacketed columns. Air-sensitive compounds have been chromatographed in an atmosphere of nitrogen or argon. All of these variations, however, are basically chromatography, and it is these basics that will be practiced in this experiment.

EXPERIMENTAL PROCEDURE

Part A

meso-Tetraphenylporphyrin, $\text{C}_{44}\text{H}_{30}\text{N}_4$ (H_2TPP)

Place a Teflon-coated stirbar in a 100-mL one-neck round-bottom flask, and add 40 mL of propanoic acid (sometimes called propionic acid). Fit the flask with a reflux condenser (see Fig. 23-5), and use a small amount of silicone grease on the joint. Leave the top of the reflux condenser open to air. Heat the acid to reflux with a rheostat-controlled heating mantle. When the propanoic acid begins to boil vigorously, add a mixture of 1.65 mL (15.75 mmol) of benzaldehyde and 1.0 mL (14.4 mmol) of pyrrole by pouring this solution down the reflux condenser. Measure out these liquids with a syringe or pipet. (*Note:* Freshly distilled pyrrole gives a higher yield of product but is not required.) Rinse the pyrrole and benzaldehyde down the condenser with 10 mL of propanoic acid. Continue to reflux the solution for 30 min, and then remove the heat and let the flask cool for a few minutes. Filter the dark brown mixture through a medium-porosity glass frit (see Figure 13-1). Rinse the mixture with a few mL of methanol until the washings are clear and purple crystals are left on the frit. Allow the crystals to dry by pulling air through them for a few minutes. Collect the purple crystals (do not scrape the




Figure 23-5
Apparatus for preparing H_2TPP .

frit too vigorously or your sample will become contaminated with powdered glass). Record the yield. Discard the wash solutions by pouring them into the containers provided for them.

Record the proton nuclear magnetic resonance (^1H NMR) spectrum of your product. To do this, place approximately 15 mg of your sample in a 5-mm diameter NMR tube and add about 1 mL of CDCl_3 (containing a small amount of tetramethylsilane, TMS, standard). The solution should fill the lower 4 cm of the NMR tube (using less CDCl_3 will lead to problems shimming the magnet). Record the spectrum from $\delta +10$ to $\delta -5$ (relative to TMS at $\delta 0$), expand and plot any complex regions, and integrate the spectrum. Now add about 3 drops of D_2O , shake the two phase CDCl_3 - D_2O mixture vigorously for about 15 seconds, and again record the NMR spectrum.

Part B

(Tetraphenylporphyrinato)copper(II), $\text{Cu}(\text{TPP})$

The reaction set up is the same as in Part A (Fig. 23-5). Place 0.1 g (0.16 mmol) of H_2TPP in a 100-mL one-neck round-bottom flask along with a stir bar. Add 20 mL of *N,N*-dimethylformamide (DMF) and stir the mixture (it is not necessary for all of the H_2TPP to dissolve). To the dark purple solution, add 0.16 g (0.80 mmol) of hydrated copper(II) acetate, $\text{Cu}(\text{O}_2\text{CCH}_3)_2 \cdot \text{H}_2\text{O}$. Fit the flask with a reflux condenser and use grease to lubricate the joint. Bring the reaction mixture to reflux with a rheostat-controlled heating mantle. Allow the reaction to proceed for 30 min.

Because the copper complex is nonfluorescent, conversion to the copper complex can be confirmed by checking for complete quenching of the porphyrin fluorescence under long-wavelength UV light. Spot some of the reaction solution on a non-fluorescing TLC plate with a Pasteur pipet, and examine the plate with a long-wavelength UV light. If the conversion is not complete (i.e., if the spot still glows red—even around the edges), add additional copper(II) acetate to the reaction solution and reflux the mixture for 10 min; then redo the spot test. Alternatively, the conversion to the copper complex can be followed by UV-vis spectroscopy. Dip a capillary tube (open at both ends) into the reaction mixture and then dip the tube containing a little of the reaction solution into a cuvette containing pure dichloromethane. Complete conversion is indicated by the disappearance of the bands at 650 and 592 nm.

When the reaction is complete, cool the reaction mixture to room temperature in an ice-water bath for 5–10 min, and then add 50 mL of distilled water to precipitate the porphyrinic material. Transfer the slurry to a separatory funnel, and extract it three times with 25 mL of dichloromethane. Collect the organic (bottom) layer each time. Discard the aqueous layer that remains in the funnel (it typically is pale blue), and then pour the combined organic extracts back into the separatory funnel. Add 50 mL of distilled water, shake the mixture, and then drain out the organic (bottom) layer. Concentrate the organic layer to dryness on a rotary evaporator. Discard the aqueous wash solution that remains in the separatory funnel.

Chromatographic Purification of Cu(TPP)

In this step, you will use TLC to determine which solvent system will allow you to separate Cu(TPP) from the H_2 TPP and chlorin impurities.

With a very small portion of the Cu(TPP) mixture, prepare a concentrated dichloromethane solution for use in the TLC trial separations. A very small amount of this solution will be spotted on the silica gel 8–10 mm from the bottom of the microscope slide. Because the best separations are achieved when this spot is between 3 and 5 mm in diameter, the dichloromethane solution should be applied with a very small capillary. These can be purchased or can be prepared by heating the middle of a melting point capillary tube over a low flame and quickly pulling the ends of the tube apart before the tube is sealed off. Scratch the capillary with a sharp file and break the tube into two applicators. Dip an applicator into the dichloromethane solution of the mixture and touch it to the TLC plate, giving a spot that is not larger than 4 mm in diameter. Allow the dichloromethane to evaporate, and then make a second application of the solution to the same spot.

Fill five screw-cap bottles (Fig. 23-3) with the trial solvents listed below to a depth of 4 mm or less, such that the level is below your spot when the TLC plate is inserted.

1. Hexane (or petroleum ether, 60–70 °C boiling fraction)
2. Toluene
3. Ethyl acetate
4. Dichloromethane
5. Acetone

Insert a spotted TLC plate into each bottle and replace the cover. Allow the solvent to rise about three-fourths of the way up the plate and then remove it and allow it to dry. Do this for each of the five TLC plates, carefully recording which solvent was used for each slide. The Cu(TPP) and the impurities are intensely colored and easily visible even at low concentrations on the plate. In other cases, however, where the compounds to be separated are weakly colored and are difficult to see, the locations of the spots can be determined either chemically (by exposing the plate to a developer such as iodine) or by examining the plate under UV light.

Make a drawing in your research notebook of the location of the spots for each of the five attempts. To establish which spot is H_2 TPP, obtain a fresh TLC plate, spotting it first with the dichloromethane solution of the mixture and then in an adjacent position with a dichloromethane solution of pure H_2 TPP. Develop the plate in one of the solvents that gave a good separation of spots, and establish which spot of the mixture is H_2 TPP by comparing it with the known H_2 TPP spot.

From the five TLC trials, you should select a solvent for the column chromatographic separation. You may select a solvent in which one of the components moves rapidly and the other more slowly. Such a solvent should give a good separation on the column. Alternatively, one might choose an initial solvent in which only one component moves while the other remains at the starting point. The first component could then be washed off the column, and then a more polar

solvent could be added to elute the other component. Regardless of your choice, it should be based on the separation achieved on the TLC plates.

Clamp a 3-cm diameter, 30-cm long chromatography column to a ring stand or other rigid support. Push a 2-cm diameter wad of glass wool to the bottom of the tube with a rod. Then pour in enough sand to give a 1-cm layer, and add enough solvent to cover the sand. In a medium sized beaker, prepare a slurry of silica gel in your selected solvent. Use enough silica gel to fill the column to a height of about 8–10 cm. Pour the slurry gently into the column and let the excess solvent drain until the meniscus is about 0.5 cm above the settled gel. (Do not allow the solvent level to drop below the top of the gel.) Carefully pour a slurry of silica gel and the Cu(TPP) mixture in a few milliliters of the initial solvent onto the column. (Depending on which solvent you choose, the sample may or may not dissolve completely. Even if it does not dissolve completely, a good separation should result if you have selected your solvent correctly.) Lower the solvent again at the top of the adsorbent and add another 1-cm layer of sand to prevent the bed from being disturbed when the eluting solvent is added. Then gently fill the column with the initial eluting solvent, being careful not to agitate the bed. Carry out the elution, using a flow rate of a few drops per second. Insoluble residues, if any, will remain on top of the column. Generally, two bands move down the column: a slowly moving sharp band that is composed of impurities, and a quickly moving broad red smear that contains the desired Cu(TPP) product. Collect this broad band continuously until all the red material has been eluted, or until the sharp band is about 2 cm from the bottom of the column. Collect only the broad red band due to Cu(TPP).

Concentrate the Cu(TPP) solution to 50 mL on a rotary evaporator and then add 50 mL of methanol to precipitate the product. Filter the solution using a medium glass frit, rinse the product with two 10-mL portions of methanol, and dry the product under vacuum if possible. Record the weight. The used silica gel can be removed by inverting the column over a large-mouth jar. The tar that remains at the top of the column should be discarded along with the used silica gel; be sure to put them into the appropriate waste container.

Obtain the visible (or UV-vis) spectrum of a very dilute toluene solution of the purified Cu(TPP) product (the molar absorption coefficient in benzene is reported to be 20,600 at 538 nm).

REPORT

Include the following:

1. Yields of H₂TPP and Cu(TPP).
2. The ¹H NMR spectra of H₂TPP both with and without added D₂O.
3. Visible (or UV-vis) spectra of H₂TPP and Cu(TPP).

PROBLEMS

1. What are the assignments of the NMR peaks for H₂TPP? Why are some of the peaks multiplets while others are singlets? Why are some of the peaks downfield of TMS (positive shift) and some upfield (negative shift)?

2. How does the NMR spectrum of H₂TPP change upon addition of D₂O? Why are only some of the peaks affected by D₂O?
3. What are the electronic transitions that are responsible for the peaks in the visible spectra of H₂TPP and Cu(TPP), and why are the two spectra different?
4. Why does Cu(TPP) elute before H₂TPP?
5. Why should the reaction of benzaldehyde and pyrrole *not* be carried out under a N₂ atmosphere?
6. In the purification of Cu(TPP), what might be the tar that remained at the top of the chromatography column?
7. The rates at which compounds elute from a silica gel column depend on any pretreatment of the silica gel. In which case would a compound elute faster: down a column made of silica gel that had previously been heated at 150 °C under vacuum for 8 h, or down a column made of silica gel that had been sitting open in the laboratory for a few days? Explain.
8. What methods might be used to detect the elution of colorless compounds from a column?
9. A mixture of *cis* and *trans* isomers of the neutral complex Cr(CO)₄[P(C₆H₅)₃]₂ is loaded onto a silica gel column and eluted with CHCl₃. Which isomer would elute first, and why?

INDEPENDENT STUDIES

- A. Prepare and characterize the iron porphyrin complex Fe(TPP)Cl and the imidazole adduct [Fe(TPP)(Im)₂]Cl, which are related to the oxidized form of the oxygen-carrying molecule hemoglobin. Use UV-vis to show isosbestic behavior in the reaction of Fe(TPP)Cl with imidazole. (Epstein, L. M.; Straub, D. K.; Maricondi, C. *Inorg. Chem.* **1967**, *6*, 1721; Collins, D. M.; Countryman, R.; Hoard, J. L. *J. Am. Chem. Soc.* **1972**, *94*, 2066.)
- B. Brominate Cu(TPP) at all eight pyrrole positions. (Bhyrappa, P.; Krishnam, V. *Inorg. Chem.* **1991**, *30*, 239.)
- C. Prepare and characterize the nickel complex Ni(TPP), which is a model of the nickel-containing enzyme called F-430. (Johnson, E. C.; Dolphin, D. *Inorg. Synth.* **1980**, *20*, 143.)
- D. Prepare the water soluble porphyrins tetrakis(carboxymethylphenyl)porphyrin or tetrakis(alkylpyridyl)porphyrin. (Datta-Gupta, N.; Bardos, T. J. *J. Heterocycl. Chem.* **1966**, *3*, 495; Beckmann, B. A.; Bochman, A.; Pasternack, R. F.; Reinprecht, J. T.; Vogel, G. C. *J. Chem. Educ.* **1976**, *53*, 387.)
- E. Investigate the use of metals as templates for the synthesis of large-ring compounds such as [Ni(R₄[14]-1,3,8,10-tetraeneN₄)]^{+ / 0}. (Tait, A. M.; Busch, D. H. *Inorg. Synth.* **1978**, *18*, 22; Hayes, J. W. II; Taylor, C. J.; Hotz, R. P. *J. Chem. Educ.* **1996**, *73*, 991.)
- F. Prepare Goedken's macrocycle and investigate its coordination chemistry. (Chipperfield, J. R.; Woodward, S. *J. Chem. Educ.* **1994**, *71*, 75.)
- G. Study the uptake of O₂ by a cobalt(II) Schiff-base complex. (Aymes, D. J.; Paris, M. R. *J. Chem. Educ.* **1989**, *66*, 854.)

- H. Use microscale techniques to prepare H_2TPP , $Zn(TPP)$, and $Ni(TPP)$. (Marsh, D. F.; Mink, L. M. *J. Chem. Educ* **1996**, *73*, 1188.)
- I. Record and interpret the EPR spectrum of $Cu(TPP)$. (For a recent EPR study of a copper(II) complex, see Louloudi, M.; Deligiannakis, Y.; Tuchagues, J.-P.; Donnadieu, B.; Hadjiliadis, N. *Inorg. Chem.* **1997**, *36*, 6335.)

REFERENCES

Porphyrins and Metalloporphyrins

- Adler, A. D.; Longo, F. R.; Finarelli, J. D.; Goldmacher, J.; Assour, J.; Korsakoff, L. *J. Org. Chem.* **1967**, *32*, 476. Synthesis of H_2TPP .
- Adler, A. D.; Longo, F. R.; Váradi, V. *Inorg. Synth.* **1976**, *16*, 213. Synthesis of H_2TPP .
- Abraham, R. J.; Jackson, A. H.; Kenner, G. W.; Warburton, D. *J. Chem. Soc.* **1963**, 853. NMR spectrum of H_2TPP .
- Collman, J. P.; Wagenknecht, P. S.; Hutchison, J. E. *Angew. Chem., Int. Ed. Engl.* **1994**, *33*, 1537. A review of face-to-face porphyrin complexes.
- Dolphin, D. *The Porphyrins*, Academic: New York, 1979. An excellent multivolume review of porphyrin chemistry.
- Lindsey, J. S.; Schreiman, I. C.; Hsu, H. C.; Kearney, P. C.; Marguerettaz, A. M. *J. Org. Chem.* **1987**, *52*, 827. Alternative synthesis of H_2TPP and other porphyrin rings.
- Scheidt, W. R.; Reed, C. A. *Chem. Rev.* **1981**, *81*, 543. Review of iron porphyrin chemistry.
- Mashiko, T.; Dolphin, D. in *Comprehensive Coordination Chemistry*, Wilkinson, G.; Gillard, R. D.; McCleverty, J. A., Eds, Pergamon: New York, 1987; Chapter 21. Excellent review of metalloporphyrins and related complexes of phthalocyanines and other macrocyclic nitrogenous ligands.

Chromatographic Techniques

- Braithwaite, A.; Smith, F. J. *Chromatographic Methods*, 4th ed., Chapman and Hall: New York, 1985.
- Gritter, R. J.; Bobbitt, J. M.; Schwarting, A. E. *Introduction to Chromatography*, 2nd ed., Holden-Day: Oakland, CA, 1985. An excellent, practical introduction to thin layer, column, and gas chromatography. Contains extensive general references and a list of chromatographic equipment suppliers.
- Druding, L. F.; Kauffman, G. B. *Coord. Chem. Rev.* **1968**, *3*, 409. Thin layer, column, and paper chromatography of coordination complexes.
- Guiochon, G.; Pommier, C. *Gas Chromatography in Inorganics and Organometallics*, Ann Arbor Science: Ann Arbor, MI, 1973. Excellent coverage of experimental techniques and the older literature.
- Miller, J. M. *Chromatography: Concepts and Contrasts*, Wiley: New York, 1988.
- Poole, C. F.; Poole, S. K. *Chromatography Today*, Elsevier: New York, 1991.

For General References to Bioinorganic Chemistry, see Experiment 21

15201). It is essential to boil the reaction solution vigorously. Reaction times up to 18 h are acceptable (and even beneficial).

Experiment 20. This experiment is fairly advanced. The $[\text{Ni}(\text{NH}_3)_6]\text{Cl}_2$ should be fresh because it loses NH_3 upon standing. This complex can be prepared while the Na reacts with the cyclopentadiene. Cyclopentadiene has an objectionable odor and is toxic, thus the cracking should be conducted in a hood. The "Cp cracker" can be reused for many weeks; any unused cyclopentadiene can be added back to the distillation flask. Do not fill the distillation flask more than about one-third full with dicyclopentadiene: The cracking process is accompanied by foaming, which can carry uncracked dicyclopentadiene into the receiver. The synthesis and purification of nickelocene requires one busy lab period followed by a lab period that is dedicated only to the sublimation step. At the end of the experiment, we often asked the students to flame-seal their samples in ampules on the vacuum line (nickelocene has a short shelf life in air). This assignment gives students some experience with glassblowing.

Experiment 21. This experiment is a useful introduction to the use of inert atmosphere techniques. Good results are obtained even by those with no previous experience with inert atmosphere techniques. We recommend that the solution of Br_2 in acetic acid be prepared in advance by the instructor.

Experiment 22. We used relatively inexpensive pH meters and standard electrodes with good results. The standardizing buffers were purchased.

Experiment 23. This experiment works well despite the low yields ($\sim 10\%$) for the H_2TPP synthesis. Pyrrole oxidizes slowly in air, and better yields of H_2TPP can be obtained if the pyrrole is freshly distilled. The experiment can be shortened by omitting the TLC study of $\text{Cu}(\text{TPP})$ and proceeding directly to the chromatographic purification using dichloromethane as the eluting solvent.

Cuvettes for UV-vis spectroscopy can be purchased from many supply houses; we used products from Hellma Cells (Box 544, Borough Hall St., Jamaica, NY 11424). A hand-held 4-watt UV light works well for detecting the free H_2TPP ; these can be purchased from Cole Parmer. For column chromatography, we used 0.040–0.063 mm, 230–400 mesh ASTM silica gel (Fisher Scientific, 1600 W. Glenlake, Itasca, IL 60142). For TLC, we used non-fluorescent Silica Gel 60 plates, which were obtained in boxes of 20×20 cm sheets with a layer thickness of 0.25 mm (EM Industries, Inc., 480 Democrat Road, Gibbstown, NJ 08027). The sheets were cut with scissors into 1×5 cm strips, which were stored in a screw-capped jar.